COLESTEROL HDL (MÉTODO DIRECTO)

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1. Gotto, A.M., Lipoprotein metabolism and the etiology of hyperlipidemia, Hospital Practice, 23; Suppl. 1, 4 (1988).
2. Crouse, J.R. et al., Studies of low density lipoprotein molecular weight in human beings with coronary artery disease, J. Lipid Res., 26; 566 (1985).
3. Badimon, J.J., Badimon, L., Fuster V., Regression of Atherosclerotic Lesions by High Density Lipoprotein Plasma Fraction in the cholesterol-Fed Rabbit, Journal of Clinical Investigation, 1990; 85: 1234-41.
4. Castelli, W.P. et al., HDL cholesterol and other lipids in coronary heart disease, Circulation, 55; 767 (1977).
5. Barr, D.P., Russ E.M., Eder, H.A., Protein-lipid relationships in human plasma, Am. J. Med., 11; 480 (1951).
6. Gordon, T. et al., High density lipoprotein as a protective factor against coronary heart disease, Am. J. Med., 62; 707 (1977).
7. Williams, P., et al., High density lipoprotein and coronary risk factor, Lancet, 1; 72, (1979).
8. Kannel, W.B., Castelli, W.P., Gordon, T., cholesterol in the prediction of atherosclerotic disease; New perspectives based on the Framingham study, Ann. Intern. Med., 90; 85, (1979).
9. National Institutes of Health publication No. 93-3095, September, (1993).
10. Special Communication, Executive Summary of the Third Report of the National cholesterol Education Program (NCEP) Expert Panel on Detection, Evaluation, and Treatment of High Blood cholesterol in Adults (Adult Treatment Panel III), JAMA, Vol. 285, No. 19, May 16, 2001, pages 2486 – 2497.
11. Warnick, G. Russell, Wood, Peter D., National Cholesterol Education Program Recommendations for Measurement of High-Density Lipoprotein Cholesterol: Executive Summary, Clinical Chemistry, Vol. 41, No. 10, 1427-1433 (1995).
12. Kimberly MM, Leary ET, Cole TG and Waymack PP. Selection, validation, standardization, and performance of a designated comparison method for HDL-cholesterol for use in the cholesterol reference method laboratory network. Clin Chem 1999; 45: 1803-12.

13. National Committee for Clinical Laboratory Standards, National Evaluation Protocols for Interference Testing, Evaluation Protocol Number 7, Vol. 6, No. 13, August (1986).
14. Young, D.S., Effects of Drugs on Clinical Laboratory Tests, 3rd ed., AACC Press, Washington, DC, 1990, 3-104 thru 3-106.
15. Carey RN, Garber CC. Evaluation of methods. In: Kaplan LA, Pesce AJ, eds. Clinical Chemistry: theory, analysis and correlation. Third Edition. St. Louis: The CV Mosby Company, St. Louis, MO., 1996:402-423.
16. Westgard JO, Carey RN, Wold S. Criteria for judging precision and accuracy in method development and evaluation. Clinical Chemistry 1974;20:825-833.