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REF 285-465

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#### **Referencias:**

- Hajna, A.A. 1955. A new specimen preservative for gram-negative organisms of the intestinal group. Public Health Lab. 13:59-62.
- Hajna, A.A. 1955. A new enrichment broth medium for gram-negative organisms of the intestinal group. Public Health Lab. 13:83-89.
- Croft, C.C., and M.J. Miller. 1956. Isolation of *Shigella* from rectal swabs with Hajna "GN" broth. Am. J. Clin. Pathol. 26:411-417.
- Taylor, W.I., and D. Schelhart. 1967. Isolation of shigellae. IV. Comparison of plating media with stools. Am. J. Clin. Pathol. 48:356-362.
- Downes, F.P. and K. Ito (ed.). 2001. Compendium of methods for the microbiological examination of foods, 4th ed. American Public Health Association, Washington, D.C.
- National Committee for Clinical Laboratory Standards. 2001. Approved Guideline M29-A2. Protection of laboratory workers from occupationally acquired infections, 2nd ed. NCCLS, Wayne, Pa.
- Garner, J.S. 1996. Hospital Infection Control Practices Advisory Committee, U.S. Department of Health and Human Services, Centers for Disease Control and Prevention. Guideline for isolation precautions in hospitals. Infect. Control Hospital Epidemiol. 17:53-80.
- U.S. Department of Health and Human Services. 1999. Biosafety in microbiological and biomedical laboratories, HHS Publication (CDC), 4th ed. U.S. Government Printing Office, Washington, D.C.
- Directive 2000/54/EC of the European Parliament and of the Council of 18 September 2000 on the protection of workers from risks related to exposure to biological agents at work (seventh individual directive within the meaning of Article 16(1) of Directive 89/391/EEC). Official Journal L262, 17/10/2000, p. 0021-0045.
- Murray, P.R., E.J. Baron, J.H. Jorgensen, M.A. Tenover, and R.H. Tenover (ed.). 2003. Manual of clinical microbiology, 8th ed. American Society for Microbiology, Washington, D.C.
- Forbes, B.A., D.F. Sahm, and A.S. Weissfeld. 2002. Bailey & Scott's diagnostic microbiology, 11th ed. Mosby, Inc., St. Louis.
- Ewing, W.H. 1986. Edwards and Ewing's identification of Enterobacteriaceae, 4th ed. Elsevier Science Publishing Co., Inc., New York.
- MacFaddin, J.F. 1985. Media for isolation-cultivation-identification-maintenance of medical bacteria, vol. 1. Williams & Wilkins, Baltimore.
- Taylor, W.I. and D. Schelhart. 1968. Isolation of *Shigellae*. V. Comparison of enrichment broths with stools. Appl. Microbiol. 16: 1383-1386.

Rev. 3: 04/2021 CIO